



Mercury Series Gas Cells
10 centimeter single-pass gas cell



Part Number

Description
Mercury Series Gas Cells

Price

MERC-0.1L/10CM-VCR-PRG

Mercury Series Single Pass Gas Cell. Baseplate and mounting fixture for center-focus FTIR with telescoping purge tubes to seal against sample compartment bulkheads (1/8" SS Swagelock purge connections provided). Sample line connections 1/4" SS welded VCR fittings. Comes with two 37.5mm x 4mm KBr windows.

\$3,800.00

Pathlength: 10 cm

Volume: 0.1 liters

Cell Body: electropolished SS

Cell Hardware: aluminum black anodized

Optical Coating: n/a


Temperature Range: ambient-200 C

Pressure Range: vacuum-30 p.s.i.

Part #: MERHJ-10CM-KIT-110

Description: Heating Kit for 10 cm Mercury Series gas cell, includes insulating hardware and heater

\$940.00



Gemini Mercury
 Series Gas Cells

